

New and improved Portable Cooler range

With our new Portable Vaccine Coolers, you have the **versatility and simplicity** to take your vaccines with you on the move.

With two sizes ranging from **8L to 18L**, you can plan your movements for the day without the worry of cold chain compromise. Alongside this, the use of gel packs rather than a refrigeration system, means you can take this product anywhere **without the worry of needing to find a plug**.



8 Litre Cooler
VCP8UK | 444410882



18 Litre Cooler
VCP18UK | 444410883



- + **+2°C to +8°C** operating temperature range
- + **External** LED temperature display
- + Holds temperature inside between 2°C & 8°C for **up to 30 hours**¹
- + Comes with **9 Gel Packs** (VPC8UK) and **12 Gel Packs** (VPC18UK)
- + Aluminium heat preserving foil lining
- + Typical vaccine storage **2.5-3L** (8L) and **12-13L** (18L)
- + Covered bung non slip feet
- + Webbing adjustable shoulder strap²

1. 8L model holds the same temperature for up to 17 hours 2. 18L model only.